

Bang Bang Shaker

UNI-EQUIP



Features:

- Digital control for speed & timer
- RPM automatically tuning function
- Auto alarm function
- Usable for incubator room
- 40 kg loading capacity
- Speed 15-350 rpm
- Orbit 45 mm

Specifications

MODEL	BBT-1	BBT-2	BB-1	BB-2
Dim. (WxDxHmm)	480 x 610 x 90	800 x 610 x 90	480 x 610 x 90	800 x 610 x 90
Installed Plate	250ml x 23Ea	250ml x 39Ea	250ml x 23Ea	250ml x 39Ea
Platform (WxD)	460 x 465	770 x 465	460 x 465	770 x 465
Control system	Time & Speed controller		Speed controller	
Speed range	15 to 350rpm			
Power	220 Volt, 60 Hz			
Weight (Kg)	32	55.3	32.5	55.8

Dimensions and specifications are subject to change without notice for further improvements on design.
Dimensions may vary dependent on different accessories be used.

UNI-EQUIP

D-82152 Martinsried / München - Fraunhoferstrasse 11
 Telefon: +49-89-8 57 52 00 • Telefax: +49-89-8 56 13 04
 e-mail: info@uniequip.de • Internet: www.uniequip.com