



# CERTOCLAV

*The Personal Sterilizer: 140°C + 2.7 bar with large capacity*



- Temperature range: 125°C (1.4 bar) & 140°C (2.7 bar)
- Made of aluminium
- 2-Step valve
- Thermometer
- Precision manometer with maximum indicator
- Safety lock (no opening under pressure)
- Steam release valve and steam release clock
- Anti-block safety valve
- Safety standard approved
- Heater included

Liter	12 l	18 l
Weight	6.1 kg	7.7 kg
Chamber diameter / Internal height	24cm / 26 cm	24 cm / 40 cm
Certoclav accessories	Instrument holder, instrument basket with tripod, condensation lid, instrument tray, tripod	

Standards which hitherto were only possible with larger sterilizers are now attainable in every suitable working situation with this practical work bench unit. The CertoClav is easy to operate achieving full sterilization within 3 minutes on the 2nd setting at 140°C (=284°F) and 2.7 bar. Setting 1 allows sterilization of heat sensitive materials within 15 minutes at 125°C (=257°F). The CertoClav is easy to transport.

With a CertoClav as their personal sterilizer qualified staff are freed from the restrictions commonly encountered when having to rely on -, and /or share, a large unit. Everything can now be sterilized at the workplace as and when required: cultures, glass containers, metal, rubber, steam-proof plastics, textiles etc. The second step of 140°C allows the full sterilization required in the bacteriological field. Often a restricted budget will not permit the purchase of a much needed sterilizing unit. The extremely competitive price of the certoclav therefore lies well within the scope of the decision making body in each laboratory.

Dimensions and specifications are subject to change without notice for further improvements on design.  
Dimensions may vary dependent on different accessories be used.



D-82152 Martinsried / München - Fraunhoferstrasse 11  
 Telefon: +49-89-8 57 52 00 · Telefax: +49-89-8 56 13 04  
 e-mail: info@uniequip.de · Internet: www.uniequip.com