

UNICHROMAT 350

UNIEQUIP

Chromatography Cabinet for small Äkta Systems

"Äkta is a trademark of GE Healthcare companies"



- Volume: 350 litres
- Interior and front panel made of stainless steel
- White housing, powder coated
- Digital display in 0,1°C
- FCKW - free cooling system, fully hermetically and with forced ventilation, mounted on anti-vibration damper, smooth running and low noise emission
- Working temperature: +4°C constant - caused on amplified cooling performance, for compensating the heating emission of the Äkta-System
- Fully insulated safety-glass doors, large glass area for best viewing conditions
- Illumination
- One EU plug-in on the inside wall left.
- blower for air circulation
- Evaporator is mounted outside of the chamber on top
- Easy cleaning
- Adjustable feet for 100% levelling
- NTC 10k resistance
- Power 220V / 50Hz

Protection of the ÄKTA System:

- Under temperature safety device with own sensor
- Over temperature protection

Optional Accessories:

- Potential-free outlet for connection to a house alarm-system
- Stainless steel shelves, adjustable
- Additional moisture-proofed plugs
- Retractable drawers
- Chromatography-column holder
- Data logger with alarm and software

Dimensions and specifications are subject to change without notice for further improvements on design.
Dimensions may vary dependent on different accessories be used.

UNIEQUIP

D-82152 Martinsried / München - Fraunhoferstrasse 11
Telefon: +49-89-8 57 52 00 • Telefax: +49-89-8 56 13 04
e-mail: info@uniequip.de • Internet: www.uniequip.com