

UNICRYO Freeze Drying Duo System

ColdT Trap with 1.360cm²/2.040 condensing surface

More than 25 years success

- Cold trap capacity 2x2/ liters
- Each trap offers 680cm² or 1.020cm² Condensing surface
- Temperature max. -55°C
- Including vacuum gauge
- Defrost heater for each condensing trap
- Round shelf racks with heat able platforms Available.
- Fuzzy logic
- CFC free
- Noiseless
- Large accessory program



The new dual UNICRYO Freeze Dryer offers a very big condensations capacity over the 2 cold traps. Each condensing trap is very fast activated by simply switch on the front panel. After occupation of ice crust and liquid evaporation of the first trap, the second trap is ready and offers new condensation of 680°C. The UNICRYO can process 2 different freeze drying jobs with its both traps. This is high tech freeze drying in compact size. UniEquip supports research with German engineering.

Dimensions and specifications are subject to change without notice for further improvements on design.
Dimensions may vary dependent on different accessories be used.