


UNICRYO

Freeze Dryer MC -60°C

UNIEQUIP

- Vessel capacity 2,4 or 6 liters
- 750 to 3500 ml freezing process per 24 hours
- Transparent drum
- Variable drying platforms with heated surface
- Vacuumgauge included
- Defrosting heater
- Small footprint
- Fuzzy logic



Dimensions and specifications are subject to change without notice for further improvements on design.
Dimensions may vary dependent on different accessories be used.

UNIEQUIP

D-82152 Martinsried / München - Fraunhoferstrasse 11
Telefon: +49-89-8 57 52 00 · Telefax: +49-89-8 56 13 04
e-mail: info@uniequip.de · Internet: www.uniequip.com