

UNICRYO COLD TRAP

25 years young with long time success

UNIEQUIP

- Cold Trap capacity
2 or 4 or 6 liters
- Condensing capacity
from 750 to 4500 ml
- Vacuumgauge included
- Max. temperature -60°C
- Defrost heater
- Fuzzy logic
- CFC free
- Noiseless



Accessories PVD Closure plate with vacuum in and outlet
GE Glass insert trap for condensing isotop material

MC 2 L -60°C	MC 4 L -60°C	MC 6 L -60°C
Danfoss 1/6 PS Compressor	Danfoss 1/5 PS	Danfoss 1/5 PS
CFC free system	CFC free system	CFC free system
Freon R507 -x 70 gr	Freon R507-x 90gr	Freon R507-x 100 gr
Power 220 / 240 V	Power 220 / 240 V or 110 V	Power 220 / 240V or 110V
Amp / Watt 1,8 / 400	Amp / Watt 2 / 450	Amp / Watt 2 / 450
Condenser surface 680 cm ²	Condenser surface 1020 cm ²	Condenser surface 1460 cm ²
Condenser made of stainless steel	Condenser made of stainless steel	Condenser made of stainless steel
Outside dim. 290 x 535 x 265	Outside dim. 400 x 685 x 255	Outside dim. 400 x 700 x 335
Weight 24 kg	Weight 32 kg	Weight 35 kg

Dimensions and specifications are subject to change without notice for further improvements on design.
Dimensions may vary dependent on different accessories be used.

UNIEQUIP

D-82152 Martinsried / München - Fraunhoferstrasse 11
 Telefon: +49-89-8 57 52 00 · Telefax: +49-89-8 56 13 04
 e-mail: info@uniequip.de · Internet: www.uniequip.com