

Laboratory Freezer -30°C

Energy class A++

UNIEQUIP

UNICRYO 180

188 liters capacity

4 drawers

Dim. h/w/d: 125/66/68.3 cm

UNICRYO 230

233 liters capacity

5 drawers

Dim. h/w/d: 144,7/66/68.3 cm

UNICRYO 280

280 liters capacity

6 drawers

Dim. h/w/d: 164,4/66/68.3 cm

UNICRYO 300

305 liters capacity

7 drawers

Dim. h/w/d: 184,1/66/68.3 cm



- Temperature range from -15°C down to -30°C
- MagicEye with digital temperature display
- Quantity-controlled automatic super frost function
- Control lamps for mains, super frost, alarm, frost control
- Visual and audible door and temperature alarm
- Interface for connection to home bus or TeleSafe module
- Frost safe system with completely enclosed, transparent drawers
- Re-circulated air cooling including Non Frost System
- Automatic defrosting & dual freeze function
- Drawers mounted on telescopic rails
- Vario space with removable drawers and intervening shelves
- Slimline door handle with integrated opening mechanism
- Flush with wall and with front ventilation
- 2 cold storage accumulators
- Adjustable feet
- Interior light

Power only in 220V / 50Hz, CE certificated, Made in GERMANY

Dimensions and specifications are subject to change without notice for further improvements on design.
Dimensions may vary dependent on different accessories be used.

UNIEQUIP

D-82152 Martinsried / München - Fraunhoferstrasse 11
Telefon: +49-89-8 57 52 00 · Telefax: +49-89-8 56 13 04
e-mail: info@uniequip.de · Internet: www.uniequip.com