

UNICRYO

Laboratory Freezer -40°C

Energy class A+++

UNIEQUIP

- Size 650 liters
- Digital display
- Alarm system
- Key locker & vent. stop
- Data acquisition system
- -40°C CFC free cooling system
- Hermetic compressor built in cascade system
- Stainless steel interior with drawers or shelf
- Outside with powder paint steel construction
- 160 mm wall isolation with special thermic splitting
- Front door with special silicon rubber>heated
- Compressor ½ PS low noise level EU made
- Drawers made in stainless steel or aluminium
- Blood Freezer with special dividers
- Data Logger temperature recorder
- 200V 50Hz



DIFFERENT DRAWERS AND INVENTORY SYSTEM AVAILABLE. Picture shows stainless steel sliding drawers for blood bags (optional).

Technical information	
Net value	~ 550liters > shelf capacity
Compressor	¾ PS > R 404 > CFC free
Automatic defrost system supported by air flow circulation	
Variable shelf levelling system on 15 positions	
Weight	135 kg net
Power	3,5 Amp > 600 Watt (5,5 kWh/24h)
External dim. WxDxT	200x76x70 cm
Internal dim. WxDxT	144x63x56 cm

Dimensions and specifications are subject to change without notice for further improvements on design.
Dimensions may vary dependent on different accessories be used.

UNIEQUIP

D-82152 Martinsried / München - Fraunhoferstrasse 11
Telefon: +49-89-8 57 52 00 · Telefax: +49-89-8 56 13 04
e-mail: info@uniequip.de · Internet: www.uniequip.com