

## UNITWIST 3-D Rocker Shaker



### For all laboratory applications

- **Low weight: 7 kg**
- **10 kg carrying capacity**
- **Shaking platform 300x300mm**
- **Moisture resistant enclosed motor**
- **Relative humidity maximum 95%**
- **Variable speed from 10 to 100 rpm**
- **Bigger shaking platforms on request**
- **Operating temperature from  $-5^{\circ}\text{C}$  up to  $+60^{\circ}\text{C}$**
- **Dimensions: 335x335x200 mm (WxDxH)**
- **Possibility for combination with UNITHERM HB water bath for hybridizations, washes, ...**

Dimensions and specifications are subject to change without notice for further improvements on design.  
Dimensions may vary dependent on different accessories be used.