

UNITWIST

Multishaker

UNIEQUIP



The UNITWIST Multishaker is supplied with a unique motion technology, which can move forwards and backwards. It is an exactly controlled cycle rotation which depends on the microprocessor controller with timer. The shaker is supplied as standard with a 260x260mm platform including silicon mat and fixing springs.

- Multi-orbital shaking motion – continuous or forward & reverse orbital shaking
- Dimensions (wxdxh) w/o platform: 260x305x100mm
- Microprocessor controlled
- Multiphase Motor (Step Down Motor)
- Settable from 0 up to 200rpm > 1 rpm increment
- Orbit 20 mm
- Max. 20kg loading cap.
- Timer 0-999 min
- Weight: 8kg

Different platforms available :

Universal platform 300x300mm for clamps
Universal platform 300x300mm with 2 spring bars, rubber coated

Universal platform 420x420mm for clamps
Universal platform 500x420mm for clamps
Universal platform with 4 spring bars, rubber coated

Available clamp sizes:
POM, autoclavable and unbreakable: 250ml, 500ml, 1000ml;
Stainless steel: 25ml, 50ml, 100ml, 250ml, 500ml, 1000ml, 2000ml;



Dimensions and specifications are subject to change without notice for further improvements on design.
Dimensions may vary dependent on different accessories be used.

UNIEQUIP

D-82152 Martinsried / München - Fraunhoferstrasse 11
Telefon: +49-89-8 57 52 00 • Telefax: +49-89-8 56 13 04
e-mail: info@uniequip.de • Internet: www.uniequip.com