


UNITWIST RT

Rocking Table Shaker



For use in cold rooms, incubators, ovens,...
Or in combination with UNIHOOD

- Compact design, smooth running, strong motor
- Dimensions: 335x335x190mm
- Adjusted tilt 5.5° (4.5° or 7.5° also available)
- Operating temperature from -5°C up to +60°C
- Platform 300x300mm with anti-slip mat
- Moisture resistant enclosed motor
- Relative humidity maximum 95%
- Variable speed from 1 to 50 rpm
- 5 kg carrying capacity
- Low weight: 7 kg

Dimensions and specifications are subject to change without notice for further improvements on design.
Dimensions may vary dependent on different accessories be used.