



UNICHROMAT

Chromatography Cabinets +4°C

- Temp. range betw. +1°C to +18°C
- Automatic defrost system
- Capacity: 900 / 1500 / 1700 liters
- Digital temperature control
- Dual Alarm system
- Lock optional
- Stainless steel chamber
- Fully isolated glass doors
- Special design on request
- We offer special stainless steel sliding drawer (Mod. 900) or ss. table for ÄKTA HPLC equipment



UNICHROMAT 900



UNICHROMAT 1500 UNIVERSAL

NEW VERSION of the 1500 Model

Since October 2015 we offer the mod 1500 w/o the middle vertical support. This version is an answer to the new generation of GE HPLC instruments with changeable instrumentation from the left to right side. Advantage: the user can work in front of a total open working area.

This type of the Unichromat is similar looking to the 1700 Model but does not offer the same depth of the chamber. Special doorstoppers fix the doors in place. Inside the cabinet we offer a stainless steel table 580 deep x 860 long x 550 mm high incl. 1 SS drawer for chemicals. On left or right side we offer movable SS Racks with 4 drawers. The temperature recovery after closing the doors back to +4°C we manage in 7 seconds

The new version of Unichromat 1500 we call
Unichromat 1500 Universal

The UNICHROMAT Series

Model

Unichromat 900

Unichromat 1500 Universal

Unichromat 1700 (AVANT)

Dimensions (W x D x H)

external dim in mm

internal dim in mm

780 x 900 x 1900

650 x 820 x 1480

1420 x 790 x 1980

1280 x 660 x 1470

1550 x 950 x 1980

1400 x 850 x 1450

Temperature Range

+1 up to +18°C
recommended

+1 up to +18°C
recommended

+1 up to +18°C
recommended



No. of Shelves

4

4

4

No. of Doors

1

2

2

Äkta HPLC equipment are manufactured by GE Healthcare - Life Science

WE SUPPLY STAINLESS STEEL DRAWERS FOR ÄKTA SYSTEMS