

Electrophoresis Chambers

Horizontal

UNI-EQUIP

Manufactured to the highest specification using the latest Injection Moulding Techniques, these units are durable, strong and economic and combine ease of use with high specification sample resolution.



MINI

| | |
|----------|----------------|
| MINI-7 | 7x7cm |
| MINI-10 | 7x10cm |
| MINI-DUO | 7x7cm & 7x10cm |

Including comb: 2x8 teeth, 1mm thickness, casting dams, loading guides.

Optional combs:

thickness: 0.75 / 1 / 1.5 / 2mm teeth: 8 / 10 / 12 / 16

MIDI

| | |
|----------|------------------|
| MIDI-7 | 10x7cm |
| MIDI-10 | 10x10cm |
| MIDI-DUO | 10x7cm & 10x10cm |

Including comb: 2x16 teeth, 1mm thickness, casting dams, loading guides.

Optional combs:

thickness: 0.75 / 1 / 1.5 / 2mm teeth: 8 / 10 / 12 / 16 / 20 / 25



CHOICE

| | |
|-------------|---------------------------|
| CHOICE-7 | 15x7cm |
| CHOICE-10 | 15x10cm |
| CHOICE-15 | 15x15cm |
| CHOICE-TRIO | 15x7cm, 15x10cm & 15x15cm |

Including comb: 2x20 teeth, 1mm thickness, casting dams, loading guides.

Optional combs:

thickness: 0.75 / 1 / 1.5 / 2mm teeth: 10 / 12 / 14 / 16 / 18 / 20 / 28 / 30 / 35



MAXI

| | |
|----------|-------------------|
| MAXI-20 | 20x20cm |
| MAXI-10 | 20x10cm |
| MAXI-DUO | 20x10cm & 20x20cm |

Including: 2x20 teeth, 1mm thickness, casting dams, loading guides.

Optional combs:

thickness: 0.75 / 1 / 1.5 / 2mm teeth: 10 / 16 / 20 / 25 / 30 / 36 / 40 / 50



Larger model: SCREEN-32 (26x32cm), SCREEN-24 (26x24cm), SCREEN-16 (26x16cm), SCREEN-TRIO (26x16cm, 24cm & 32cm).

Dimensions and specifications are subject to change without notice for further improvements on design.
Dimensions may vary dependent on different accessories be used.

UNI-EQUIP

D-82152 Martinsried / München - Fraunhoferstrasse 11
Telefon: +49-89-8 57 52 00 · Telefax: +49-89-8 56 13 04
e-mail: info@uniequip.de · Internet: www.uniequip.com