



UNI-EQUIP

Electrophoresis Chambers

Vertical

PAGE MINI

Plate dimensions 10x10cm
Gel dimensions 7.5x8cm

Max. sample capacity 80 samples / 20 samples per gel

Buffer volume Min. 250ml / Max. 1200ml

Including: 2 set of glass plates, 1mm thick bonded spacers, 2x12 teeth, 1mm thick combs.

Combs available:

No. of teeth 1, 5, 8, 9, 10, 12, 16, 20
Thickness 0.75, 1, 1.5, 2mm

A wide range of accessories and different inserts are available !

This Mini Vertical unit is injection moulded making it one of the most economic Mini Vertical units currently available. Easy casting and gel loading ensure rapid, trouble free set up times.

PAGE MAXI

Plate dimensions 20x20cm
Gel dimensions 16x17.5cm

Max. sample capacity 192 samples / 48 samples per gel

Buffer volume Min. 1200ml / Max. 5600ml

Including: 2 set of glass plates, 1mm thick bonded spacers, 2x12 teeth, 1mm thick combs.

Combs available :

No. of teeth 1, 5, 10, 18, 24, 30, 36, 48
Thickness 0.75, 1, 1.5, 2mm

A wide range of accessories and different inserts are available !



Casting base

This Maxi vertical unit for large vertical gel separations is based on the same design and enhances resolution in extended separations.

We supply DAIICHISILVERSTAIN II Kits !

Dimensions and specifications are subject to change without notice for further improvements on design.
Dimensions may vary dependent on different accessories be used.

UNI-EQUIP

D-82152 Martinsried / München - Fraunhoferstrasse 11
Telefon: +49-89-8 57 52 00 · Telefax: +49-89-8 56 13 04
e-mail: info@uniequip.de · Internet: www.uniequip.com