

# Portable pH-Meter



Model PHB-4 (ATC)

Model PHB-1 (MTC)

## Model PHB-1 Portable pH Meter

- Intelligent chip design
- Automatic calibration and manual temperature compensation

Model	PHB-1
Measuring range	0.00 to 14.00pH +/-1999mV
Resolution	0.01pH 1mV
Accuracy	+/-0.01pH +/-%full scale
Temperature comp.	0 to 80°C (MTC)
Electrode	201-A plastic pH combination electrode
Buffer solution	PH 4.00, 7.00, 10.01/btl.
Power	DC9V battery
Size and weight	168x43x26mm/170g

## Model PHB-4 Portable pH/mV/Temp Meter

- Intelligent chip design
- Automatic calibration and automatic temperature compensation

Model	PHB-4
Measuring range	0.00 to 14.00pH +/-1999mV 0 to 80°C
Resolution	0.01pH 1mV 0.1°C
Accuracy	+/-0.01pH +/-%full scale +/-1°C
Temperature comp.	0 to 80°C (ATC)
Electrode	201B-S pH/ATC trinal combination electrode
Buffer solution	PH 4.00, 7.00, 10.01/btl.
Power	DC9V battery
Size and weight	168x43x26mm/170g