

INCUBATOR TO GO ANYWHERE...



UNI-CUB MINI INCUBATOR

New range of portable, top loading mini incubators

UNI-CUB Mini Incubator is a compact, easy to use incubator designed to operate consistently, even in the toughest conditions. The distinctive "top loading" design allows quick, easy access to the chamber without having to disturb the contents.

The unit is manufactured from high quality plastics which are not only tough and durable but lightweight for easy handling; the UNI-CUB Mini 10 weighs 1800 grammes. The DC powered UNI-CUB Mini 20 can be powered from a car cigarette lighter plus the mains for maximum flexibility.



The easy to use controls make it ideal for use by non technical staff. The only operator controls are an on/off switch and a numbered dial for setting the temperature. The raised lid includes four ventilation plugs. These can be removed to allow increased air circulation to avoid a build up of condensation within the chamber.

All UNI-CUB Mini instruments, will accurately hold the set temperature within +/- 1°C and heat up from 20°C to 35°C in only 20-30 minutes.



	UNI-CUB 10	UNI-CUB 20
Voltage	230 VAC (115 VAC available to order)	12 VDC (Inc Power Supply 90-264 VAC)
Frequency	50-60 Hz	DC
Temperature Range	Ambient +5°C to 50°C	Ambient +5°C to 47°C
Temperature Accuracy	Mean ±1.0°C at single point	
Temperature Variation	Mean ±2.0°C within Chamber	
Power Rating	80 watts	45 watts
Fuse Rating	T315mA	F5A
Weight	1800 grams	2700 grams
Gross Weight	2.0 kgs	3.0 kgs
Controls	Main on/off switch, Heating Controller & Indicator Fuse Holder on rear panel	
Dimensions	Internal Chamber: 180 x 150 x 175(H) mm Overall External: 320 x 250 x 260(H) mm Max Internal Volume: 5 litres	
Packed Dimensions	380 x 310 x 300 mm	
Capacity	The incubator's chamber will comfortably hold 20 Dipslides, 20 Petri Dishes or 50 Contact Plates.	

Applications

Due to the compact yet tough design, the UNI-CUB Mini units are ideally suited to use on site, out of the laboratory, and particularly with Dipslides. The prime areas of use for this instrument are in the industrial, environmental and food safety sectors.

UNI-CUB Mini incubators can be supplied with either a Raised Clear Lid or a Raised White Lid. Alternatives are available on request.

A location is provided within the chamber for the use of a thermometer.

All units are CE marked and supplied fitted with plugs, a manual & a 12 Month return to base warranty.

See you @ uniequip.com...

