

Thermal Cycler

Peltier-based Thermal Cycler

MG 16+

MG 25+



- **High Quality and Long Life**
- **Heated Lid**
- **Small, Light and Noiseless**
- **LED display**
- **Peltier driven**

Technology	Peltier-base
Block Capacity	0.2ml×25well/0.5ml×16well
Block Temperature Range	0°C ~ 99°C ,but not more than 30°C below ambient temp.
Block Heating Rate	3°C/s
Block Cooling Rate	2°C/s
Block Uniformity	< ±0.5°C
Block Accuracy	< ±0.5°C
Mode of Temperature Control	Sim-Tube & Block
No. of Programs	62
Max. No. of Cycles	99
Temp. Increment/Decrement	1 ~ 60s
Time Increment/Decrement	1 ~ 10°C
Auto Restart	Yes
Hold at 4°C	Forever
Running Time Display	Yes
Hot Cover Temperature	105°C
Dimension(L×W×H)	325×220×135mm
Net Weight	3.6kg
Wide Power Supply Range	110~220V 50/60Hz 150W