

MICROCENTRIFUGE UEC-Force 6 *Speedy Egg*



Features

- Personal benchtop microcentrifuge
- Angle rotor is applicable for:
 - 6x 1.5/ 0.5/0.2 tubes
 - Or two 8x 0,2ml strip tubes
- Easy operation close lid, push start and it runs
- Open lid - the rotor stops immediately or press stop
- Extraordinarily stable while running

Specifications	UEC-Force 6
Max. Speed	6000 rpm
Max. RCF	2000 xg
Weight	~ 800 g
Dimensions	189x161x121mm
Rotor (inc.)	6x 1.5 / 0.5 / 0.2ml tubes Two 8x 0,2ml strip tubes or 16x individual 0,2ml tubes

www.uniequip.com

Dimensions and specifications are subject to change without notice for further improvements on design.
Dimensions may vary dependent on different accessories be used.

UNI-EQUIP

D-82152 Martinsried / München - Fraunhoferstrasse 11
Telefon: +49-89-8 57 52 00 · Telefax: +49-89-8 56 13 04
e-mail: info@uniequip.de · Internet: www.uniequip.com

High-Speed Mini Centrifuge UEC 14



Features

- Low noise level – extremely quiet
- Maintenance free drive
- LCD display – easy to follow
- Capacity: 12 x 1.5/2.0ml tubes
- Separate short spin key
- Sample heating only 12°C after 20 min. at max. speed
- Automatic lid release
- Standard rotor made of anodized aluminium
- Autoclavable rotor (121°C, 20min.)
- Max. speed of 14.500 rpm
- Max. RCF of 14.000 x g
- RMP/RCF setting as required
- Timer can be set to 99 min.

Model	UEC 14
Power supply	100-240V ~ 50-60Hz
Power	105W
Max. speed	14.500 rpm
Max. centrifugal force	14.000 x g
Timer	20 sec. to 99 min.
Max. load	12 x 1.5/2.0ml tubes
Max. permitted density	1.2g/ml
Acceleration time to max. speed	20 sec.
Braking time from max. speed	19 sec.
Ambient temperature	5°C – 35°C
Dim. WxDxH	212mm x 240mm x 138mm
Weight (incl. rotor)	4.8kg

www.uniequip.com

Dimensions and specifications are subject to change without notice for further improvements on design.
Dimensions may vary dependent on different accessories be used.

UNI-EQUIP

D-82152 Martinsried / München - Fraunhoferstrasse 11
 Telefon: +49-89-8 57 52 00 · Telefax: +49-89-8 56 13 04
 e-mail: info@uniequip.de · Internet: www.uniequip.com