

UNIMAG - Talent

Magnetic stirring + heating

UNIEQUIP

This instrument series combines magnetic stirring up to a range of 1800rpm (adjustable) with a 18x18cm heating plate to reach up to 450°C*+ > real time temperature adjustment by means of an external temperature probe. Tripod included.

- Microprocessor controller
- Large LCD screen to display more values at same time
- Surface specially coated for chemical resistance
- External temperature probe for medium
- Free-step rotation for speed adjustment
- Excellent suited for viscous liquids
- Over temperature safety device
- Heat indication for caution
- Platform size 18x18cm



Model	Stirring in Liter capacity)	Speed range (rpm)	Liquid temp range (°C)	Working plate temp. (°C)	Temp. display resolution	Plate dimension (mm)	Consumption (W)
IT-UE-09A5	5	100-1800	RT~ 250	--	±1%	180x180	420
IT-UE-09A8	8	100-1800					630
IT-UE-09A12	12	100-1800					1040
IT-UE-09B5	5	100-1800	--	RT~320	±1%		420
IT-UE-09B10	10	100-1800		RT~380			630
IT-UE-09B15	15	100-1800		RT~450			1040
External dimension: 210x330x130mm / Weight between 3,5 & 4,5 kg							

* depends on the model

Dimensions and specifications are subject to change without notice for further improvements on design.
Dimensions may vary dependent on different accessories be used.

UNIEQUIP

D-82152 Martinsried / München - Fraunhoferstrasse 11
 Telefon: +49-89-8 57 52 00 · Telefax: +49-89-8 56 13 04
 e-mail: info@uniequip.de · Internet: www.uniequip.com