

UNITWIST 3-D

Rocker Shaker



For all laboratory applications

- **Low weight: 7 kg**
- **10 kg carrying capacity**
- **Shaking platform 300x300mm**
- **Moisture resistant enclosed motor**
- **Relative humidity maximum 95%**
- **Variable speed from 10 to 100 rpm**
- **Bigger shaking platforms on request**
- **Operating temperature from -5°C up to $+60^{\circ}\text{C}$**
- **Dimensions: 335x335x200 mm (WxDxH)**
- **Possibility for combination with UNITHERM HB water bath for hybridizations, washes, ...**

UNITWIST RT

Rocking Table Shaker



For use in cold rooms, incubators, ovens,...
Or in combination with UNIHOOD

- **Compact design, smooth running, strong motor**
- **Dimensions: 335x335x190mm**
- **Adjusted tilt 5.5° (4.5° or 7.5° also available)**
- **Operating temperature from -5°C up to +60°C**
- **Platform 300x300mm with anti-slip mat**
- **Moisture resistant enclosed motor**
- **Relative humidity maximum 95%**
- **Variable speed from 1 to 50 rpm**
- **5 kg carrying capacity**
- **Low weight: 7 kg**