

Unicryo

Ice Flake Machine



- Air cooled compressor works free of FCKW (on request water cooled)
- Inlet for operation with distilled water optionally available
- Housing made in stainless steel, storage container in plastic
- Electronic control

Model	ES50A	ES70A	ES100A	ES200A
Ice output per 24h	50 kg	72 kg	100 kg	200 kg
Ice storage	12 kg	17 kg	25 kg	47 kg
Dimensions cm	420x600x720	500x600x810	630x640x960	930x560x1090
Height incl. feet	+3,5cm	+8cm	+8cm	+8cm