

CERTOCLAV

*The Personal Sterilizer: 140°C + 2.7 bar
with large capacity*



- Temperature range: 125°C (1.4 bar) & 140°C (2.7 bar)
- Made of **aluminium**
- 2-Step valve
- Thermometer
- Precision manometer with maximum indicator
- Safety lock (no opening under pressure)
- Steam release valve and steam release **clock**
- Anti-block safety valve
- Safety standard approved
- Heater included

Liter	12 l	18 l
Weight	6.1 kg	7.7 kg
Chamber diameter / Internal height	24cm / 26 cm	24 cm / 40 cm
Certoclav accessories	Instrument holder, instrument basket with tripod , condensation lid, instrument tray , tripod	

Standards which hitherto were only possible with larger sterilizers are now attainable in every suitable working situation with this practical work bench unit. **The CertoClav** is easy to operate achieving full sterilization within 3 minutes on the 2nd setting at 140°C (=284°F) and 2.7 bar. Setting 1 allows sterilization of heat sensitive materials within 15 minutes at 125°C (=257°F). **The CertoClav** is easy to transport.

With a CertoClav as their personal sterilizer qualified staff are freed from the restrictions commonly encountered when having to rely on -, and /or share, a large unit. Everything can now be sterilized at the workplace as and when required: cultures, glass containers, metal, rubber, steam-proof plastics, textiles etc. The second step of 140°C allows the full sterilization required in the bacteriological field. Often a restricted budget will not permit the purchase of a much needed sterilizing unit. The extremely competitive price of the certoclav therefore lies well within the scope of the decision making body in each laboratory.