

SA-300 H 40 liters chamber

Automatic Autoclav Sterilizer



- **Wide range of usage**
This autoclav provides sterilization temperature from 121°C – 132°C. It is **suitable for** wrapped / unwrapped materials, utensils, liquids, and etc. Widely used in hospital, clinic and laboratories.
- **Low Water Alarm**
When the water in the chamber or in the water reservoir is not enough, the “low water” indication lamp will light on, buzzer and stop all functions at the same time.
- **Pre-Heating-Program**
Pre-heating by dry heater, to keep the chamber temperature in blance. This leads to the perfect dry function and best sterilization

Model	SA – 300 H
Chamber size	Stainless steel #304 300mm Ø x 550mm (D)
Overall	600(W) x 470(H) x 740(D)mm
Power consumption	224V / 2.5KW / 10.5A 220V / 2.1KW / 9.5A
Using temperature / Pressure	Unwrapped: 121-132°C / 2.1kg / cm ² Wrapped: 121-132°C / 2.1kg / cm ² Liquids: 121°C / 1.2kg / cm ²
Design temperature	140°C
Safety devices	Pressure control switch, door safety catcher, over heat protection, safety valve, electronic circuit safety system, low water indication and cut-off, emergency exhaust switch
Standard accessories	Stainless steel #304 Sterilization box X 1 pc. Tray X 1 pc.
Safety approval	CE – EMC, LVD, by TÜV CE-MDD approval by DNV Chamber designed by ASME/USA standard
Auto-Dry Function	Yes
Chamber capacity	40 liters