

Bang Bang Shaker



Features:

- Digital control for speed & timer
- RPM automatically tuning function
- Auto alarm function
- Usable for incubator room
- 40 kg loading capacity
- Speed 15-350 rpm
- Orbit 45 mm

Specifications

MODEL	BBT-1	BBT-2	BB-1	BB-2
Dim. (WxDxHmm)	480 x 610 x 90	800 x 610 x 90	480 x 610 x 90	800 x 610 x90
Installed Plate	250ml x 23Ea	250ml x 39Ea	250ml x 23Ea	250ml x 39Ea
Platform (WxD)	460 x 465	770 x 465	460 x 465	770 x 465
Control system	Time & Speed controller		Speed controller	
Speed range	15 to 350rpm			
Power	220 Volt, 60 Hz			
Weight (Kg)	32	55.3	32.5	55.8



Shaker

It gains the rotational force using both poles of electro magnetic force. It is an ultra power saving type shakes that added powerful rotational force by installing the fly wheel of metal cast mold.

Optional plates are for following sizes available >>>

100ml x 36, 100ml x 72, 250ml x 20, 250ml x 39,
500ml x 16, 500ml x 32, 1000ml x 9, 1000ml x 16

Racks of multi Shaker:

Model MBB-1,
Model MBB-2,

Dim. (WxDxH): 680 x 640 x 1950mm
Dim. (WxDxH): 1000 x 640 x 1950mm