


UNITHERM WAM 6/15T

Perspex-waterbath
incl. hang-in-thermostate
+ magnetic stirrer



Waterbath

Allround dimensions: 596x270x187mm
Bath dimensions: 410x254x80mm
Weight: 6,2kg

Hang-in-Thermostate

Working temp. 25°C to 100°C
Temp. permanence: +/- 0,3k
Temp. setting: digital
Forwarding amount: 5l / Min.
Heating power: 1,0 KW
Safety class: Class II– DIN12879
Weight: 3,3 kg

Magnetic-Stirrer UNIMAG 15/6

Adjustable from: 50 to 850l / Min.

Max. stirring capacity:
UNIMAG 15 15x 250ml
UNIMAG 6 6x 500ml

Performance:
UNIMAG 15 15 W
UNIMAG 6 8 W

Weight:
UNIMAG 15 2,8kg
UNIMAG 6 2,2kg

