

Laboratory Freezer -30°C

Energy class A++

UNICRYO 180
188 liters capacity
4 drawers
Dim. h/w/d: 125/66/68.3 cm

UNICRYO 230
233 liters capacity
5 drawers
Dim. h/w/d: 144,7/66/68.3 cm

UNICRYO 280
280 liters capacity
6 drawers
Dim. h/w/d: 164,4/66/68.3 cm

UNICRYO 300
305 liters capacity
7 drawers
Dim. h/w/d: 184,1/66/68.3 cm



- Temperature range from -15°C down to -30°C
- MagicEye with digital temperature display
- Quantity-controlled automatic super frost function
- Control lamps for mains, super frost, alarm, frost control
- Visual and audible door and temperature alarm
- Interface for connection to home bus or TeleSafe module
- Frost safe system with completely enclosed, transparent drawers
- Re-circulated air cooling including Non Frost System
- Automatic defrosting & dual freeze function
- Drawers mounted on telescopic rails
- Vario space with removable drawers and intervening shelves
- Slimline door handle with integrated opening mechanism
- Flush with wall and with front ventilation
- 2 cold storage accumulators
- Adjustable feet
- Interior light

Power only in 220V / 50Hz, CE certificated, Made in GERMANY