

UEP-Pipette


The new generation micropipette

- High performance, light weight and easy in easy operating
- Fatigue free working with convenient finger rest
- Made with special material for better protection against chemical and physical corrosion
- Pipette handle has coating of Thermoplastic Elastomer to prevent transfer of body heat to internal parts for high accuracy even on continuous use
- Volume display is always visible while pipetting
- In house cleaning, repair and calibration possible without any special tools
- New tip cone design for smooth fitting with most of manufacturer's tips
- Color coding for easy identification
- Fully autoclavable without disassembly for maximum convenience as well as protection from contamination
- Enhanced UV resistance
- Shelf mounting stand for efficient handling and storage
- Supplied with individual quality certificate (QC) and calibration report according to ISO8655
- CE IVD Compliant



- Manifold rotates 360° to allow pipetting from any convenient angle
- Durable V-Rings for smooth loading and ejection of variety of tips with minimum force as well as for leak free operation
- The tip ejector is designed with convenient shape to reduce tip loading and ejection forces
- Includes all basic features of the single pipette





- Light & smooth plunger action
- Locking mechanism preventing accidental volume changes during pipetting



- Clear 4 digit display for small increments and wider volume selection options
- Resistance free click stop counter
- Larger digits



- PTFE sealing is used to create precise air tightness to provide long term consistent performance



- Slim design allows pipette to operate even with narrow tubes used in laboratory
- Easily accessible recalibration mechanism without any chance of accidental change in calibration



- Ergonomic designing
- Suitable for left and right handed users with a relaxed grip and a good balance
- Soft touch tip ejection



UEP Variable Volume Micropipette			UEP Fix Volume Micropipette		UEP 8 Channel Micropipette		
Cat. No.	Volume [µl]	Increm. [µl]	Cat. No.	Volume [µl]	Cat. No.	Volume [µl]	Increm. [µl]
UEP-V-125	0,1-0,5	0,002	UEP-F-205	2,5	UEP-8-050	5-50µl	0,100
UEP-V-100	0,5-10	0,020	UEP-F-5	5	UEP-8-100	10-100µl	0,100
UEP-V-200	5-50	0,100	UEP-F-10	10	UEP-8-200	20-200µl	0,100
UEP-V-300	50-200	0,200	UEP-F-20	20	UEP-8-300	30-300µl	0,200
UEP-V-400	200-1000	1,000	UEP-F-25	25	UEP-8-010	0,5-10µl	0,020
UEP-V-500	10-100	0,100	UEP-F-50	50	UEP 12 Channel Micropipette		
UEP-V-600	100-1000	1,000	UEP-F-100	100	Cat. No.	Volume [µl]	Increm. [µl]
UEP-V-700	2-20	0,020	UEP-F-200	200	UEP-12-050	5-50µl	0,100
UEP-V-800	20-200	0,200	UEP-F-250	250	UEP-12-100	10-100µl	0,100
UEP-V-900	1000-5000	10,000	UEP-F-500	500	UEP-12-200	20-200µl	0,100
UEP-V-905	500-5000	10,000	UEP-F-1000	1000	UEP-12-300	30-300µl	0,200
UEP-V-1000	1000-10000	20,000	UEP-F-2000	2000	UEP-12-010	0,5-10µl	0,020
UEP-V-1005	500-10000	20,000	UEP-F-5000	5000	UEP Pipette Stand		
			UEP-F-10000	10000	Cat. No.	Description	
					UEP-ST-01	Carousel stand for 6 pipettes	
					UEP-ST-02	Shelf mounting stand for individual pipette	