

Horizontal Table Top Autoclave

CLASS B – TÜV and CE certificated, ISO 9001



- Pre-Vacuum & Dry-Vacuum
- CLASS B – TÜV and CE certificated, ISO 9001
- 16 or 24 liters chamber capacity
- Sterilization temperature 121°C / 135°C
- Intelligent design for your professional choice
- Installed 24 programs and special programs for PRION sterilization, LIQUID sterilization and test programs, full automatic processing
- **Auto-fill** water system
- Door lock indication
- Digital temperature display
- Analog vacuum gauge

Model	SA-260MA	SA-230MA
Auto dry time can be selected. 0, 15, 30, 35, 40 mins	Yes	Yes
Ext dimensions (dxwxh)	655mm x 533mm x 442mm	594mm x 522mm x 410mm
Water capacity of water tank	4200cc	4200cc
Chamber capacity	24 liters	16 liters
Net weight	52 kg	49 kg
Power supply	230V 50/60Hz	230V 50/60Hz
Power consumption	2400W / 11Amp.	1900W / 8Amp.
SAFETY DEVICE		
Pressure safety valve & protection switch	Yes	Yes
Over temperature protection	Yes	Yes
Water tank water level sensor	Yes	Yes
Pressure door auto-lock	Yes	Yes
Standard accessories	Sterilization Box, Tray 1 pc.	Sterilization Box, Tray 1 pc.
OPTIONAL ACCESSORIES	Sterilization Plate set	Sterilization Plate set