

UNIGELDRYER

Model 3545



DRYING YOUR GELS AS FAST AS POSSIBLE

- Digital temperature and time control
- Aluminium made drying surface with dimensions: 35x45 cm
- Scintered polyethylene sheet
- Sealed drying area by silicone cover
- New heat cover and vacuum lid
- Effective and quick heat transfer
- Fast drying by pre-vacuum
- 100% corrosion resistant
- Steel housing
- SODA-active-prefilter optional
- Possible in RACK-SYSTEM with vacuum manifold fits for 3 UNIGELDRYERS and 1 vacuum pump