

UNIMAG ZX³



- Vortex mixers for mixing contents of test tubes or small flasks.
- The ZX³ operates on "manual" or "continuous" mode.
- Structure made in aluminium casting, epoxy painted in order to give to the instrument an excellent resistance to the attack of chemical and mechanical agents and to corrosion in general.
- Very strong sucker feet for strong fixing and high stability of the apparatus and to avoid vibrations transmission to the working bench.

Pick-up attachment:	Ø55 mm
Stirring frequency:	0-40 Hz
Power:	45 W
Weight:	2 kg

Accessories:

- Mounting for Eppendorf tubes 25x 1,5/2,0ml
- Mounting for microtitre plate
- Supporting plate Ø90 mm