

# UNITWIST 200

Small footprint – High performance



The new redesigned UNITWIST 200 is made for virtually any vortexing and mixing application. It is equipped with a basic platform with compability to all microtitre plates. Each platform contains 2 microtitreplates. Contemporary 4 plates usable. With a maximum shaking speed of 1200 upm, a justable orbit between 3 – 4mm and an integrated vortexing funktion, the UNITWIST 200 ensures rapid and complete mixing of even the smallest volumes. Its compatible with fast changeable micro test tubes, PCR plates, deepwell plates and MTPs. Speed is variable and is easily to setup. Shaking times can be set between 0-30 minutes, continuously run is possible. The UNITWIST 200 with its humidity resistance up to 85% is ideal for placement inside an cold storage room or incubator. It features an operating concept that is as simple as intelligent.

## Specifications:

- for microplates / 2 plates per platform
- 2 platforms mountable
- for 56 reaction vials 1,5/2 ml
- for microtubes 2x48 reaction vials 0,2/0,4/0,5 µl
- basic platform with:
  - anti-slip mat
  - with springbars
  - with clamps

Uniequip supports research.