

VORTEMP

Shaking Incubator

UNIEQUIP



- Programmable speed, temperature and time
- Ramps for temperature and time settable in 5 steps
- Large digital display
- Temperature range from 37°C up to 85°C (optional temp. up to 100°C)
- Large capacity: 56 microtubes or 4 microplates
- Usable also in cold rooms
- Easy interchangeable platforms for microtubes and microplates
- RS232 for GMP control
- Optional control software for PC
- Small footprint

Dimensions and specifications are subject to change without notice for further improvements on design.
 Dimensions may vary dependent on different accessories be used.

UNIEQUIP

D-82152 Martinsried / München - Fraunhoferstrasse 11
 Telefon: +49-89-8 57 52 00 · Telefax: +49-89-8 56 13 04
 e-mail: info@uniequip.de · Internet: www.uniequip.com



UNIEQUIP

Features of the new control:

- **2 Modes: Normal operation and programmable operation**
- **5 Steps programmable**
- **128x24 Dot Matrix Large Screen**
- **Soft-Start timer**
- **Speed settable in rpm**

Vortemp offers maximum of precision, speed, reliability and high capacity at minimum space consumption.

Specifications:

Temp. range:	RT +5°C up to +80°C (optional up to +100°C)
Temp. control:	microprocessor controlled
Orbit:	3mm
Speed:	250 – 1350rpm
Capacity:	56 tubes 1,5/2,0ml, 192 tubes 0,5µl, 0,2ml tubes (60pcs.), four 96 well microtiter plates ;
Heating power:	300 Watt
Dimensions (WxDxH):	280x240x210mm
Weight:	9 kg

Accessories:

- **Platform for 56 tubes 1,5/2,0ml**
- **Platform for 192 tubes 0,5µl**
- **Adapter for 0,2ml tubes (60 pcs.)**
- **Platforms (two each) for four 96 well microtiter plates**

Dimensions and specifications are subject to change without notice for further improvements on design.
Dimensions may vary dependent on different accessories be used.

UNIEQUIP

D-82152 Martinsried / München - Fraunhoferstrasse 11
Telefon: +49-89-8 57 52 00 · Telefax: +49-89-8 56 13 04
e-mail: info@uniequip.de · Internet: www.uniequip.com